

## Schisantherin A

Catalog No: #S2913


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Schisantherin A
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	536.57
Formulation	C30H32O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C[C@H]1CC2=CC3=C(C(=C2C4=C(C(=C(C=C4[C@@H]1(C)O)OC(=O)C5=CC=CC=C5)OC)OC)OC)OCOC3Pathway  
s:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zhou E, et al. Int Immunopharmacol. 2014 Sep;22(1):133-40.

Note: This product is for in vitro research use only