

CDDO(Bardoxolone)

Catalog No: #S2915


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CDDO(Bardoxolone)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Nrf2 inhibitor
Calculated MW	491.7
Formulation	C31H41NO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@@]12CC[C@]3(CCC(C[C@H]3[C@H]1C(=O)C=C4[C@]2(CC[C@@H]5[C@@]4(C=C(C(=O)C5(C)C)C#N)C)(C)C)C(=O)O

Pathways:Immunology/InflammationReceptor:Nrf2Boling pt:Melting pt:Solubility:Soluble in DMSOAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Reisman SA, et al. J Am Soc Nephrol. 2012 Oct;23(10):1663-73.

Note: This product is for in vitro research use only