

## Patchouli alcohol


Catalog No: #S2916

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Patchouli alcohol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	222.2
Formulation	C15H26O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC1=C[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC=C4[C@@]3(CC[C@@H](C4)O)C)C  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Solid  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yu JL, et al. J Surg Res. 2015 Apr;194(2):537-43.

Note: This product is for in vitro research use only