

## Cholesteryl Acetate

Catalog No: #S2918


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cholesteryl Acetate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	428.7
Formulation	C29H48O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:C[C@H](CCCC(C)C)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC=C4[C@@]3(CC[C@@H](C4)OC(=O)C)C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zanoun A, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2005 Nov;62(1-3):547-51.

Note: This product is for in vitro research use only