

## Sesamin

Catalog No: #S2919


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sesamin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	354.36
Formulation	C <sub>20</sub> H <sub>18</sub> O <sub>6</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: C1[C@H]2[C@H](CO[C@@H]2C3=CC4=C(C=C3)OCO4)[C@H](O1)C5=CC6=C(C=C5)OCO6 Pathways: Others Receptor: Others Binding pt: Melting pt: Solubility: DMSO: 70 mg/mL warmed (197.54 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Ma L, et al. Biochem Biophys Res Commun. 2015 May 29;461(2):230-6.

Note: This product is for in vitro research use only