

Berbamine dihydrochloride

Catalog No: #S2920


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Berbamine dihydrochloride
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	680.24
Formulation	C37H42Cl2N2O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: CN1CCC2=CC(=C3C=C2C1CC4=CC=C(C=C4)OC5=C(C=CC(=C5)CC6C7=C(O3)C(=C(C=C7CCN6C)OC)OC)O)OC
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Jin X, Wu Y. *Anat Rec (Hoboken)*. 2014 May;297(5):802-9.

Note: This product is for in vitro research use only