

Sinomenine hydrochloride


Catalog No: #S2921

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sinomenine hydrochloride
Brief Description	Inhibitors
Purification	99.79%
Target Name	NF-κB inhibitor
Calculated MW	365.85
Formulation	C19H24ClNO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CN1CC[C@@]23CC(=O)C(=C[C@@H]2[C@@H]1CC4=C3C(=C(C=C4)OC)O)OC Pathways: NF-κB Receptor: NF-κB
Inhibitor
g pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Song L, et al. Biochem Biophys Res Commun. 2015 Aug 28;464(3):705-10.

Note: This product is for in vitro research use only