

Schizandrin A

Catalog No: #S2926


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Schizandrin A
Brief Description	Inhibitors
Purification	99.49%
Target Name	Others
Calculated MW	416.51
Formulation	C24H32O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@@H]1CC2=CC(=C(C(=C2C3=C(C(=C(C=C3[C@@H]1C)OC)OC)OC)OC)OC)OCPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White Sheet SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang CP, et al. J Physiol Biochem. 2014 Sep;70(3):735-47.

Note: This product is for in vitro research use only