

Trigonelline hydrochloride

Catalog No: #S2927


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trigonelline hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	173.6
Formulation	C7H8ClNO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: Cl.C[N+](C)=CC=CC(=O)C([O-])=O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhang DF, et al. Asian Pac J Trop Med. 2015 Aug;8(8):651-4.

Note: This product is for in vitro research use only