

L-Rhamnose monohydrate

Catalog No: #S2928


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L-Rhamnose monohydrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	182.17
Formulation	C6H14O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O.CC1OC(O)C(O)C(O)C1OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: > 10mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pardo E, Orejas M. Microb Cell Fact. 2014 Nov 22;13:161.

Note: This product is for in vitro research use only