

Ethyl Vanillate


Catalog No: #S2929

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ethyl Vanillate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	196.2
Formulation	C10H12O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCOC(=O)C1=CC(OC)=C(O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Namazi MR Md, Shotorbani AK Md. Iran J Med Sci. 2015 Nov;40(6):478-84.

Note: This product is for in vitro research use only