

Alphalipoic acid


Catalog No: #S2930

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alphalipoic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	206.32
Formulation	C ₈ H ₁₄ O ₂ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(=O)CCCC1CCSS1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Tayebati SK, et al. Clin Exp Hypertens. 2015 Aug 19:1-9.

Note: This product is for in vitro research use only