

## Ginsenoside Rh1

Catalog No: #S2932


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ginsenoside Rh1
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	638.87
Formulation	C36H62O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:CC(=CCC[C@@](C)([C@H]1CC[C@@]2([C@@H]1[C@@H](C[C@H]3[C@]2(C[C@@H]([C@@H]4[C@@]3(CC[C@@H](C4(C)C)O)C)O)[C@H]5[C@@H]([C@H]([C@H]([C@@H]([C@H](O5)CO)O)O)C)O)C)O)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Luo SL, et al. Lett Appl Microbiol. 2015 Jan;60(1):72-8.

Note: This product is for in vitro research use only