

Bilirubin

Catalog No: #S2934


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bilirubin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	584.68
Formulation	C33H36N4O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CC1=C(NC(=C1CCC(=O)O)CC2=C(C(=C(N2)/C=C\3/C(=C(=O)N3)C=C)C)CCC(=O)O)/C=C\4/C(=C(C(=O)N4)C=C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: SolubleAppearance:Red PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wagner KH, et al. Clin Sci (Lond). 2015 Jul;129(1):1-25.

Note: This product is for in vitro research use only