

Phenethylamine

Catalog No: #S2935


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phenethylamine
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	121.18
Formulation	C8H11N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: NCCC1=CC=CC=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Le Roux G, et al. Drug Alcohol Depend. 2015 Sep 1; 154:46-53.

Note: This product is for in vitro research use only