

# 3,4,5-Trimethoxyphenol

Catalog No: #S2936

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	3,4,5-Trimethoxyphenol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	184.19
Formulation	C9H12O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Product Description

Research Area: antioxidant SMILES: COC1=CC(O)=CC(OC)=C1OC Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Parmar VS, et al. Acta Chem Scand. 1996 Jun;50(6):558-60.

Note: This product is for in vitro research use only