

## 4'-Demethylepipodophyllotoxin

Catalog No: #S2937


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	4'-Demethylepipodophyllotoxin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Topoisomerase inhibitor
Calculated MW	400.38
Formulation	C21H20O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: COC1=CC(=CC(=C1O)OC)[C@H]2[C@@H]3[C@H](COC3=O)[C@@H](C4=CC=C(C=C24)OCO5)O Pathways: DNA Damage/DNA Repair Receptor: Topo II Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Xiao L, et al. Eur J Med Chem. 2014 Jun 10;80:267-77.

Note: This product is for in vitro research use only