

Piceid

Catalog No: #S2938


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Piceid
Brief Description	Inhibitors
Purification	99.50%
Target Name	Phospholipase inhibitor
Calculated MW	390.39
Formulation	C ₂₀ H ₂₂ O ₈
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: C₁=CC(=CC=C₁/C=C/C₂=CC(=CC(=C₂)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O₃)CO)O)O)O)OP

Pathways: Metabolism Receptor: PLA₂ Boiling pt: Melting pt: Solubility: DMSO: 78 mg/mL (199.8 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Shiyu S, et al. BMC Cell Biol, 2011, 25, 12:31.

Note: This product is for in vitro research use only