

Forskolin

Catalog No: #S2939


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Forskolin
Brief Description	Inhibitors
Purification	99.25%
Target Name	AChR antagonist
Calculated MW	410.51
Formulation	C22H34O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

SystemSMILES:CC(=O)O[C@H]1[C@H]([C@@H]2[C@]([C@H](CCC2(C)C)O)([C@@]3([C@@]1(O[C@@](CC3=O)(C)C=C)O)C)O)Pathways:NeuroscienceReceptor:Adenylyl cyclase (AC)Boiling pt:Melting pt:Solubility:Ethanol: 37 mg/mL (90.13 mM); DMSO: 82 mg/mL (199.75 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Insel PA, et al. Cell Mol Neurobiol, 2003, 23(3), 305-314.

Note: This product is for in vitro research use only