

Psoralen

Catalog No: #S2942


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Psoralen
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	186.17
Formulation	C11H6O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O=C1OC2=CC3=C(C=CO3)C=C2C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Yang L, Sun X, Geng X. Pak J Pharm Sci. 2015 Mar;28(2 Suppl):667-70.

Note: This product is for in vitro research use only