

Tanshinone IIA-sulfonic sodium

Catalog No: #S2946


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tanshinone IIA-sulfonic sodium
Brief Description	Inhibitors
Purification	99.04%
Target Name	Others
Calculated MW	396.39
Formulation	C ₁₉ H ₁₇ O ₆ SNa
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: [Na+].CC1=C(OC2=C1C(=O)C(=O)C1=C2C=CC2=C1CCCC2(C)C)S([O-])(=O)=O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Red Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Weng YS, et al. Am J Chin Med. 2015 Nov 30:1550089.

Note: This product is for in vitro research use only