

Ethyl vanillin acetate

Catalog No: #S2948


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ethyl vanillin acetate
Brief Description	Inhibitors
Purification	99.73%
Target Name	Others
Calculated MW	208.21
Formulation	C ₁₁ H ₁₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCOC1=C(OC(C)=O)C=CC(C=O)=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Gray Solid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Seebaluck R, et al. J Ethnopharmacol. 2015 Jan 15;159:137-57.

Note: This product is for in vitro research use only