

Schisandrin

Catalog No: #S2950


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Schisandrin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	432.51
Formulation	C24H32O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:COC1=C(OC)C(OC)=C2C(CC(C)C(C)(O)CC3=CC(OC)=C(OC)C(OC)=C23)=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Guo CX, et al. Pharmazie. 2015 Jan;70(1):29-32.

Note: This product is for in vitro research use only