

Camphor


Catalog No: #S2952

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Camphor
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	152.24
Formulation	C10H16O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC1(C)C2CCC1(C)C(=O)C2Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Willetts A, Kelly DR. Microbiology. 2014 Aug;160(Pt 8):1783-94.

Note: This product is for in vitro research use only