

Ligustrazine hydrochloride

Catalog No: #S2953


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ligustrazine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	172.66
Formulation	C ₈ H ₁₂ N ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: Cl.CC1=NC(C)=C(C)N=C1
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Gao HJ, et al. Neural Regen Res. 2015 May;10(5):832-40.

Note: This product is for in vitro research use only