

## Gibberellic acid

Catalog No: #S2955


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Gibberellic acid
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	346.38
Formulation	C19H22O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:C[C@@]12[C@H](C=C[C@@]3([C@@H]1[C@@H]([C@]45[C@H]3CC[C@](C4)(C(=C)C5)O)C(=O)O)OC2=O)OPathways:OthersReceptor:steroidogenesisBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Premalatha R, et al. Andrologia. 2014 Oct;46(8):902-9.

Note: This product is for in vitro research use only