

Cedryl acetate

Catalog No: #S2956


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Cedryl acetate |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 264.41 |
| Formulation | C17H28O2 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:OtherSMILES:CC1CCC23CC1C(C)(C)C2CCC3(C)OC(C)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Transparent CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sultan S, et al. Eur J Med Chem. 2013 Apr;62:764-70.

Note: This product is for in vitro research use only