

Arctigenin

Catalog No: #S2957


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Arctigenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	372.41
Formulation	C21H24O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: COC1=C(C=C(C=C1)C[C@H]2COC(=O)[C@@H]2CC3=CC(=C(C=C3)O)OC)OC Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wu X, et al. *Biochem Pharmacol.* 2015 Aug 15;96(4):323-36.

Note: This product is for in vitro research use only