

Deacetyltaxol

Catalog No: #S2960


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Deacetyltaxol |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 811.87 |
| Formulation | C45H49NO13 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research

Area: OtherSMILES: CC1=C2[C@H](C(=O)[C@@]3([C@H](C[C@@H]4[C@@]([C@H]3[C@@H]([C@@]([C2(C)C)(C[C@@H]1OC(=O)[C@@H]([C@H](C5=CC=CC=C5)NC(=O)C6=CC=CC=C6)O)O)OC(=O)C7=CC=CC=C7)(CO4)OC(=O)C)O)C)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

- Grobosch T, et al. J Anal Toxicol. 2012 Jan-Feb;36(1):36-43.

Note: This product is for in vitro research use only