

Notoginsenoside R1

Catalog No: #S2961


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Notoginsenoside R1
Brief Description	Inhibitors
Purification	98%+
Target Name	ERK inhibitor
Calculated MW	933.13
Formulation	C47H80O18
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CC(=CCC[C@@](C)([C@H]1CC[C@@]2([C@@H]1[C@@H](C[C@H]3[C@]2(C[C@@H]([C@@H]4[C@@]3(CC[C@@H](C4(C)C)O)C)O[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)O[C@H]6[C@@H]([C@H]([C@@H](CO6)O)O)C)O)O[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)CO)O)O)O)CPATHWAYS:MAPKReceptor:ERK1; ERK2Boiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- Xu Y, et al. Am J Chin Med. 2014;42(4):799-816.

Note: This product is for in vitro research use only