

Dihydrocholesterol

Catalog No: #S2962


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dihydrocholesterol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	388.68
Formulation	C27H48O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@H](CCCC(C)C)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC[C@@H]4[C@@]3(CC[C@@H](C4)O)C)CPathways:Other
sReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang X, et al. Atherosclerosis. 2015 Sep;242(1):77-86.

Note: This product is for in vitro research use only