

Cholic Acid

Catalog No: #S2963


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cholic Acid
Brief Description	Inhibitors
Purification	99.92%
Target Name	Others
Calculated MW	408.58
Formulation	C ₂₄ H ₄₀ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

systemSMILES: C[C@H](CCC(=O)O)[C@H]1CC[C@@H]2[C@@]1([C@H](C[C@H]3[C@H]2[C@@H](C[C@H]4[C@@]3(CC[C@H](C4)O)C)O)C

Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wu Y, et al. J Biomed Nanotechnol. 2015 Jul;11(7):1247-60.

Note: This product is for in vitro research use only