

Astragalus polyphenols

Catalog No: #S2964


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Astragalus polyphenols
Brief Description	Inhibitors
Purification	99.77%
Target Name	Others
Calculated MW	406.38
Formulation	C20H22O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C1=CC(=CC=C1/C=C/C2=CC(=CC(=C2O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Faint Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ma ZJ, et al. Zhongguo Zhong Yao Za Zhi. 2015 Jun;40(12):2325-9.

Note: This product is for in vitro research use only