

Deoxycholic acid

Catalog No: #S2965


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Deoxycholic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	392.57
Formulation	C ₂₄ H ₄₀ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: C[C@H](CCC(=O)O)[C@H]1CC[C@@H]2[C@@]1([C@H](C[C@H]3[C@H]2CC[C@H]4[C@@]3(CC[C@H](C4)O)C)O)C
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Walker P, Lee D. J Cosmet Dermatol. 2015 Mar;14(1):33-9.

Note: This product is for in vitro research use only