

Stigmasterol


Catalog No: #S2967

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Stigmasterol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	412.69
Formulation	C29H48O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CC[C@H](/C=C/[C@@H](C)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC=C4[C@@]3(CC[C@@H](C4)O)C)C(C)CPat
hways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please
warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kim YS, et al. BMB Rep. 2014 Aug;47(8):433-8.

Note: This product is for in vitro research use only