

## Lawsone

Catalog No: #S2969


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Lawsone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	174.16
Formulation	C10H6O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OC1=CC(=O)C(=O)C2=CC=CC=C12Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Munt DJ, et al. AAPS PharmSciTech. 2015 Dec;16(6):1425-33.

Note: This product is for in vitro research use only