

## Phytic acid


Catalog No: #S2971

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phytic acid
Brief Description	Inhibitors
Purification	1.1mol/l
Target Name	Others
Calculated MW	660.04
Formulation	C6H18O24P6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:OP(O)(=O)OC1C(OP(O)(O)=O)C(OP(O)(O)=O)C(OP(O)(O)=O)C(OP(O)(O)=O)C1OP(O)(O)=OPathways:OthersReceptor:Others  
 Boling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Lv Y, et al. Neurosci Lett. 2015 Jun 15;597:132-6.

Note: This product is for in vitro research use only