

o-Veratraldehyde


Catalog No: #S2975

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	o-Veratraldehyde
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	166.18
Formulation	C9H10O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=CC=CC(C=O)=C1OC
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Peralta MA, et al. Acta Crystallogr C. 2007 Jan;63(Pt 1):o68-72.

Note: This product is for in vitro research use only