

L(-)-sorbose

Catalog No: #S2976


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L(-)-sorbose
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	180.16
Formulation	C6H12O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C([C@@H]([C@H]([C@@H](C(=O)CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Mandlaa, et al. J Ind Microbiol Biotechnol. 2015 Jun;42(6):897-904.

Note: This product is for in vitro research use only