

## Acetylvainillin

Catalog No: #S2977


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Acetylvainillin
Brief Description	Inhibitors
Purification	99.45%
Target Name	Others
Calculated MW	194.18
Formulation	C10H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:COC1=C(OC(C)=O)C=CC(C=O)=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Faint Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Shakeel F, et al. Food Chem. 2015 Aug 1;180:244-8.

Note: This product is for in vitro research use only