

Mogroside V

Catalog No: #S2978


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mogroside V
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	1287.43
Formulation	C60H102O29
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C[C@H](CC[C@H](C(C)(C)O)O[C@H]1[C@@H]([C@H]([C@@H]([C@H](O1)CO[C@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)O)O)[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)[C@H]4CC[C@@]5([C@@]4[C@H]([C@@]6([C@H]5CC=C7[C@H]6CC[C@@H](C7(C)C)O[C@H]8[C@@H]([C@H]([C@@H]([C@H](O8)CO[C@H]9[C@@H]([C@H]([C@@H]([C@H](O9)CO)O)O)O)O)C)O)C)CP

athways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Xu F, et al. J Pharm Biomed Anal. 2015 Nov 10;115:418-30.

Note: This product is for in vitro research use only