

DL-Menthol


Catalog No: #S2979

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DL-Menthol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	156.15
Formulation	C10H20O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C[C@@H]1CC[C@H]([C@@H](C1)O)C(C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang DL, et al. J Agric Food Chem. 2002 Jan 16;50(2):262-5.

Note: This product is for in vitro research use only