

## Chondroitine sulfate

Catalog No: #S2980


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Chondroitine sulfate
Brief Description	Inhibitors
Purification	99.87%
Target Name	Others
Calculated MW	463.37
Formulation	C13H21NO15S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

SystemSMILES:CC(=O)N[C@@H]1[C@H]([C@H]([C@H](O[C@H]1O)OS(=O)(=O)O)O)[C@H]2[C@@H]([C@H]([C@@H]([C@H](O2)C(=O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Kakizaki I, et al. Carbohydr Polym. 2015 May 5;121:362-71.

Note: This product is for in vitro research use only