

Xanthone

Catalog No: #S2982


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Xanthone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	196.21
Formulation	C13H8O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O=C1C2=CC=CC=C2OC2=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 1 mg/mL warmed (5.09 mM); DMSO: 39 mg/mL warmed (198.77 mM)Appearance:Dusty Pink PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Itamany EE, et al. Bioorg Med Chem Lett. 2014 Nov 1;24(21):4939-42.

Note: This product is for in vitro research use only