

## 3,4,5-Trimethoxycinnamic acid

Catalog No: #S2983


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	3,4,5-Trimethoxycinnamic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	238.24
Formulation	C12H14O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: COC1=CC(=CC(=C1OC)OC)/C=C/C(=O)O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Lee Cl, et al. Arch Pharm Res. 2013 Oct;36(10):1244-51.

Note: This product is for in vitro research use only