

Scopoletin


Catalog No: #S2984

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scopoletin
Brief Description	Inhibitors
Purification	99%+
Target Name	AChR antagonist
Calculated MW	192.17
Formulation	C10H8O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: COC1=C(O)C=C2OC(=O)C=CC2=C1
Pathways: Neuroscience
Receptor: mAChR
Boiling pt: Melting pt:
Solubility: DMSO: 10 mM
Appearance: White Needle Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lee HI, Lee MK. Toxicol Lett. 2015 Sep 17;237(3):210-8.

Note: This product is for in vitro research use only