

## Epigallocatechin

Catalog No: #S2987


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Epigallocatechin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	306.27
Formulation	C15H14O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:C1[C@H]([C@H](OC2=CC(=CC(=C21)O)O)C3=CC(=C(C(=C3)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Wang N, et al. J Agric Food Chem. 2015 Jan 26

Note: This product is for in vitro research use only