

4-Hydroxybenzoic acid


Catalog No: #S2989

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Hydroxybenzoic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	138.12
Formulation	C7H6O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(=O)C1=CC=C(O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Farhoosh R, et al. Food Chem. 2016 Mar 1;194:128-34.

Note: This product is for in vitro research use only