

Phenibut

Catalog No: #S2990


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phenibut
Brief Description	Inhibitors
Purification	99.50%
Target Name	GABA Receptor inhibitor
Calculated MW	215.68
Formulation	C10H13NO2HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system SMILES: Cl.NCC(CC(O)=O)C1=CC=CC=C1 Pathways: Membrane transporter/Ion channel;

Neuroscience Receptor: GABAR Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Bagmetova VV, et al. Bull Exp Biol Med. 2015 May;159(1):48-52.

Note: This product is for in vitro research use only