

3,3'-Diindolylmethane


Catalog No: #S2991

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3,3'-Diindolylmethane
Brief Description	Inhibitors
Purification	99.50%
Target Name	Androgen Receptor antagonist
Calculated MW	246.31
Formulation	C17H14N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
SMILES: C(C1=CNC2=CC=CC=C1)C1=CNC2=CC=CC=C1
Pathways: Endocrinology/Hormones
Receptor: Androgen Receptor
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Prabodh K, et al. Mol Pharmacol 78:297-309, 2010

Note: This product is for in vitro research use only